



REF ASK-04020-MIH1

4_45 catheter length

inch dia. spring-wire

Radial Artery Catheterization Kit

Contents:

- Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm)
- Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- Syringe: 3 mL Luer-Lock
- Foil Pack 70% Alcohol Prep Pad
- 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Frepp® Solution
- Towel: 13" x 18" (33 cm x 46 cm)
 Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- Needle Holder
- Scissors
- Medication Label: 1% Lidocaine

- Gauze Pad: 2" x 2" (5 cm x 5 cm)
 Gauze Pad: 4" x 4" (10 cm x 10 cm)
 Dressing: OpSite®3 2-3/4" x 2-1/4" (7 cm x 6 cm)
- Dressing: BIOPATCH®4
- Suture: 3-0 Silk with Curved Needle
- ¹A trademark of Becton, Dickinson and Company.

- ²A registered trademark of CareFusion or one of its subsidiaries. ³A registered trademark of T.J. Smith & Nephew, Limited.
- ⁴A registered trademark of Johnson & Johnson Corporation

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 - 25°C (68 - 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



LOT



























REF ASK-04020-MIH1

4.45 cm catheter length

Spring-wire guide

