



REF AK-09800

035 spring-wire guide 8.5^{Fr.} 10^{cm}_{sheath}length

Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/
- Side Port and Tissue Dillator Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J" Tip on 1:
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) Injection Needle: 25 Ga. x 1" (2.54 cm) Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW 1:
- 1: 1:
- Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Syringe: 3 mL Luer-Lock 1:
- 1:
- 2:
- Syringe: 5 mL Luer-Slip 5 mL 1% Lidocaine HCI Solution and Alcohol Prep (non-sterile solution)
- 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®1 One-Step Solution with Hi-Lite Orange™ Tint 1: (non-sterile solution)
- 1: 1:
- Obturator: 8 Fr. Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile
 SharpsAway® Disposal Cup
 Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
 Drape: 24" x 36" (61 cm x 91 cm)
 Scalpel: #11
 Gauze Pad: 2" x 2" (5 cm x 5 cm)
 Gauze Pad: 4" x 4" (10 cm x 10 cm)
 Suture: 3-0 Braided Silk with Straight Needle

- ¹A registered trademark of CareFusion or one of its subsidiaries

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

This product contains dry natural rubber.

Store between 20 – 25°C (68 – 77°F). Avoid freezing and excessive heat above 40°C (104°F). Sterilized using ethylene oxide unless otherwise indicated in the contents list. Fluid path components are non-pyrogenic.



Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA





LOT

Ω



35 spring-wire guide



Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters