





1 0 cm sheath

.035 inch dia. spring-wire guide

Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

Contents:

- Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/
 Side Port and Tissue Dilator
 Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J"
 Tip on Other) with Arrow Advancer
 Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle

- Carneter: 18 Ga. x 2-1/2* (6.35 cm) Radiopaque over 20 Ga. RW introducer Needle Pressure Transduction Probe Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe

- Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe Obturator: 8 Fr.
 Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter Dust Cap: Non-Vented SharpsAway® Disposal Cup
 Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
 Drape: 24" x 36" (61 cm x 91 cm)
 Scalpel: #11
 Stopcock: 4-way High-Flow
 Gauze Pad: 2" x 2" (5 cm x 5 cm)
 Gauze Pad: 4" x 4" (10 cm x 10 cm)
 Suture: 3-0 Braided Silk with Straight Needle

- 5: 1:
- Suture: 3-0 Braided Silk with Straight Needle

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

This Product Contains Dry Natural Rubber.

Fluid path components are non-pyrogenic



LOT

























10 cm sheath length

35 inch dia. spring-wire guide

