PB-10002-000



REF MTO-09903-TS (IPNXXXXX)



## Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters

Contents:

Sheath: 9 Fr. x 3-3/8" (8.5 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/ Side Port and Tissue Dilator Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) Injection Needle: 25 Ga. x 12 (2.54 cm) Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW Syringe: 3 mL Luer-Lock Syringe: 5 mL Luer-Slip 1:

- Syringe: 5 mL Luer-Slip Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL Obturator: 8 Fr. Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter Dust Cap: Non-Vented

- Dust Cap: Non-Vented SharpsAway® Disposal Cup Drape: 21\* x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive Scalpel: #11 Stopcock: 4-way High-Flow Prep Sponge Swab Gauze Pad: 4" x 4" (10 cm x 10 cm)

- 5:
- Suture: 3-0 Braided Silk with Straight Needle Tray: Prep 1: 1:

## Rx only.

## Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic





LOT LotNumber YYYY-MM-DD





Arrow International, Inc. Subsidiary of Teleflex Incorporated

**)feleflex®** 

2400 Bernville Road Reading, PA 19605 USA

> www.teleflex.com LBL015909 (01/20)





9<sup>Fr.</sup> 8.5<sup>cm</sup><sub>sheath</sub><sub>length</sub>



Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters