



REF CS-22854-E

8.5^{Fr.} 4 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Arrow-Howes™ Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Arrowg+ard Blue® Catheter

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

- fr** Jeu de cathétérisme veineux central à quatre lumières Arrow-Howes avec cathéter Blue FlexTip Arrowg+ard Blue
- de** Arrow-Howes-Besteck für vierlumigen, zentralen Venenkatheter mit Blue-FlexTip-Arrowg+ard-Blue-Katheter
- it** Set per cateterismo della vena centrale a quattro lumi Arrow-Howes con catetere Blue FlexTip Arrowg+ard Blue
- pt** Conjunto de Cateterização Venosa Central Arrow-Howes com Lúmen Quádruplo e Cateter Blue FlexTip, Arrowg+ard Blue
- es** Conjunto de cateterización venosa central de cuatro luces Arrow-Howes con catéter Blue FlexTip Arrowg+ard Blue
- sv** Arrow-Howes set för frykanalig central venkatetrering med Blue FlexTip Arrowg+ard Blue-kateter

¹Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.


California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.39	3200
Medial 1 (14 Ga.)	0.47	5400
Medial 2 (18 Ga.)	0.36	1600
Proximal (18 Ga.)	0.38	1800

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL027972 (02/20)



REF CS-22854-E 8.5^{Fr.} 4 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Arrow-Howes™ Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Arrowg+ard Blue® Catheter