



REF SA-14703

7^{Fr.} **3**^{Lumen}

30 cm catheter length

.032 inch dia. spring-wire guide

Three-Lumen CVC

Contents:

- Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 30 cm
- Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End "J" Tip on Other) with Arrow Advancer with ECG mark
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW

- Introducer Needle: Ecnogenic 18 Ga. x 2-1/2" (6.35 cm) X I W
 Pressure Transduction Probe
 Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
 Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
 SecondSite ™ Adjustable Hub: Fastener
 SecondSite ™ Adjustable Hub: Catheter Clamp
 Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with
 adhesive adhesive
 Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive
- Safety Scalpel: #11
- Prep Sponge Swab

 Valve: MicroCLAVE™¹ Neutral Displacement Connector

¹A trademark of ICU Medical, Inc.

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (16 Ga.)	0.54	2025	
Medial (18 Ga.)	0.45	896	
Proximal (18 Ga.)	0.48	955	

^{*} Priming volumes are approximate and are done without accessories.



LOT



M















Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL050073 R01 (2023-09)









.032 inch dia. spring-wire guide



[†] Flow rate values are approximate and are determined using deionized water at 100 cm head height.