

REF SA-25853

8.5^{Fr.} 3^{Lumen}



Arrowg+ard Blue[®] Three-Lumen CVC


Contents:

- Arrowg+ard Blue® Three-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 20 cm 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with
- ECG Mark 1:
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm

- SecondSite™ Adjustable Hub: Fastener SecondSite™ Adjustable Hub: Catheter Clamp
- Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with 1:
- adhesive Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive 2:
- 1: Safety Scalpel: #11
- Prep Sponge Swab Valve: MicroCLAVE™1 Neutral Displacement Connector 2: 3:
- ¹A trademark of ICU Medical, Inc.
- All components are CE 2797 unless otherwise noted.


Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



ARROW

* Priming volumes are approximate and are done without accessories. + Flow rate values are approximate and are determined using deionized water at 100 cm head height.







8.5^{Fr.} 3^{Lumen}





Arrowg+ard Blue[®] Three-Lumen CVC