



REF **AK-25123-CDC** (IPNXXXXXX)

12 Fr. **3** Lumen **20** cm catheter length **.035** inch dia. spring-wire guide

Arrowgard Blue® Large-Bore Three-Lumen CVC for High Volume Infusions

Contents:

- 1: Arrowgard Blue® Three-Lumen Catheter: 12 Fr. (4.1 mm OD) x 20 cm origin:xxxxxxxxxxxxxxxxxxxx
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer origin:xxxxxxxxxxxxxxxxxxxx
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle origin:xxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Pressure Transduction Probe origin:xxxxxxxxxxxxxxxxxxxx
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe origin:xxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe origin:xxxxxxxxxxxxxxxxxxxx
- 1: Syringe: 3 mL Luer-Lock origin:xxxxxxxxxxxxxxxxxxxx
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm origin:xxxxxxxxxxxxxxxxxxxx
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep origin:xxxxxxxxxxxxxxxxxxxx
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®¹ One-Step Solution with Hi-Lite Orange™ Tint origin:xxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® II Locking Disposal Cup origin:xxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® Disposal Cup origin:xxxxxxxxxxxxxxxxxxxx
- 1: Drape: 41" x 55" with 4" fenestration origin:xxxxxxxxxxxxxxxxxxxx
- 1: Needle Holder origin:xxxxxxxxxxxxxxxxxxxx
- 1: Safety Scalpel: #11 origin:xxxxxxxxxxxxxxxxxxxx
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Suture: 3-0 Silk with Curved Needle origin:xxxxxxxxxxxxxxxxxxxx

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65
⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Contraindications: The Arrowgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.40	2457
Medial (12 Ga.)	1.3	10151
Proximal (12 Ga.)	1.2	11039

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Set Product of xxx
Packaged in xxx

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com
LBL051034 R02 (2024-05)



REF **AK-25123-CDC**

12 Fr. **3** Lumen **20** cm catheter length **.035** inch dia. spring-wire guide



Arrowgard Blue® Large-Bore Three-Lumen CVC for High Volume Infusions