



REF **AK-24703-ASK** (IPNXXXXXX)

**7** Fr. **3** Lumen **30** cm catheter length **.032** inch dia. spring-wire guide

# Arrowgard Blue® Three-Lumen CVC

**Contents:**

- 1: Arrowgard Blue® Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 30 cm origin:xxxxxxxxxxxxxxxxxxxx
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer origin:xxxxxxxxxxxxxxxxxxxx
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle origin:xxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Pressure Transduction Probe origin:xxxxxxxxxxxxxxxxxxxx
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe origin:xxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe origin:xxxxxxxxxxxxxxxxxxxx
- 1: Syringe: 3 mL Luer-Lock origin:xxxxxxxxxxxxxxxxxxxx
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm origin:xxxxxxxxxxxxxxxxxxxx
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep origin:xxxxxxxxxxxxxxxxxxxx
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>1</sup> One-Step Solution with Hi-Lite Orange™ Tint origin:xxxxxxxxxxxxxxxxxxxx
- 3: Dust Cap: Non-Vented origin:xxxxxxxxxxxxxxxxxxxx
- 1: SecondSite™ Adjustable Hub: Fastener origin:xxxxxxxxxxxxxxxxxxxx
- 1: SecondSite™ Adjustable Hub: Catheter Clamp origin:xxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® II Locking Disposal Cup origin:xxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® Disposal Cup origin:xxxxxxxxxxxxxxxxxxxx
- 1: Drape: 41" x 50" with 4" fenestration, with adhesive origin:xxxxxxxxxxxxxxxxxxxx
- 1: Safety Scalpel: #11 origin:xxxxxxxxxxxxxxxxxxxx
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Suture: 3-0 Braided Silk with Straight Needle origin:xxxxxxxxxxxxxxxxxxxx

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**

⚠ **WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

**Contraindications: The Arrowgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.**



Lumen

Priming Volume\* (mL)

Gravity Flow Rate† (mL/hr)

Distal (16 Ga.)	0.6	2025
Medial (18 Ga.)	0.5	921
Proximal (18 Ga.)	0.5	1004

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)LotNumber

Set Product of xxx  
Packaged in xxx

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex**  
www.teleflex.com

LBL051089 R02 (2024-04)



REF **AK-24703-ASK**

**7** Fr. **3** Lumen

**30** cm catheter length

**.032** inch dia. spring-wire guide

# Arrowgard Blue® Three-Lumen CVC

