

1032 inch dia. spring-wire guide

Arrowg+ard Blue® Two-Lumen CVC

- Arrowg+ard Blue® Two-Lumen Catheter: 8 Fr. (2.8 mm OD) x 20 cm Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End "J" Tip on Other) with Arrow Advancer with ECG Mark
- Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW Syringe: 5 mL Luer-Slip

- Sylinge. 3 mix Luer-Sip Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm SecondSite™ Adjustable Hub: Fastener SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 2797 unless otherwise noted.


Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.









Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (14 Ga.)	0.8	5237	
Proximal (14 Ga.)	0.8	3511	

Priming volumes are approximate and are done without accessories.



Representative and Importer:



Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL051236 R01 (2023-10)







REF EU-25802-N











[†] Flow rate values are approximate and are determined using deionized water at 100 cm head height.