



## **REF SA-16702**

7<sup>Fr.</sup> 2<sup>Lumen</sup>


\_032 inch dia. spring-wire guide 20<sup>cm</sup><sub>catheter</sub> length

## **Two-Lumen CVC**

## Contents:

- Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm 1:
- Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End "J" Tip on Other) with Arrow Advancer with ECG Mark
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle 1:
- Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW 1:
- 1:
- 1: 1:
- 1:
- Introducer Needle: Echogenic 18 Ga. x2-1/2" (6.35 cm) XTW Pressure Transduction Probe Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm SecondSite ™ Adjustable Hub: Fastener SecondSite ™ Adjustable Hub: Catheter Clamp Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with adhesive 1: 1:
- adhesive Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive 2:
- 1: Safety Scalpel: #11
- Prep Sponge Swab Valve: MicroCLAVE™1 Neutral Displacement Connector 2: 2:
- <sup>1</sup>A trademark of ICU Medical, Inc.
- All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.




\* Priming volumes are approximate and are done without accessories.
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



**REF SA-16702 Two-Lumen CVC** 

 $\sim$ 

 $\Box$ 



EU Authorized Representative and Importer

ECREP Teleflex Medical IDA Business and Technology Park Dublin Road, Athlone 2797 Co. Westmeath, Ireland

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

**Teleflex**®

LBL051257 R01 (2023-09)

ARROW



