



REF SA-16553

3 cm catheter length

Three-Lumen CVC

- Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 13 cm
- Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J" Tip on Other) with Arrow Advancer
- Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW

- Syringe: 5 mL Luer-Slip
 Syringe: 5 mL Luer-Slip
 Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
 Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
 SecondSite™ Adjustable Hub: Fastener
 SecondSite™ Adjustable Hub: Catheter Clamp
- Safety Scalpel: #11
 Valve: MicroCLAVETM1 Neutral Displacement Connector

¹A trademark of ICU Medical, Inc.

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



| Lumen | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) | |
|-------------------|----------------------------|-------------------------------------|--|
| Distal (20 Ga.) | 0.32 | 988 | |
| Medial (22 Ga.) | 0.38 | 370 | |
| Proximal (22 Ga.) | 0.34 | 422 | |

^{*} Priming volumes are approximate and are done without accessories.





Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL051290 R01 (2023-09)









13 cm catheter length • 018 inch dia. spring-wire guide





[†] Flow rate values are approximate and are determined using deionized water at 100 cm head height.