

- 1.

- 1:

- Injection Needle: SafetyGlide<sup>™2</sup> 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
  Syringe: 3 mL Luer-Lock
  Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
  5 mL 1% Lidocaine HCI Solution and Alcohol Prep
  3 mL Applicator 2% CHG and 70% IPA ChloraPrep®3 One-Step Solution with Hi-Lite Orange™ Tint
  SecondSite<sup>™</sup> Adjustable Hub: Fastener
  SecondSite<sup>™</sup> Adjustable Hub: Catheter Clamp
  SharpsAway® II Locking Disposal Cup

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. California Prop. 65 California Prop. 65 & WARNING: Cancer and Reproductive Harm. www.P65Warnings.ca.gov Not made with natural rubber latex.

Store between 20 - 25°C (68 - 77°F). Sterilized using ethylene oxide unless otherwise indicated in the contents

list

Fluid path components are non-pyrogenic Contraindications: The Arrowg+ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (16 Ga.)	0.38	2490	
Medial (18 Ga.)	0.37	1162	
Proximal (18 Ga.)	0.39	1328	

\* Priming volumes are approximate and are done without accessories.
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

Â	MD	i			STERING	STERILE EO	)	Ť	(	LAREX	20°C (68°F)	- 25°C (77°F)
---	----	---	--	--	---------	------------	---	---	---	-------	----------------	------------------

LOT

 $\nabla$ 

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

> **Teleflex**® LBL051348 R01 (2022-02)

[لمم

**REF AK-42703-AX** 

7<sup>Fr.</sup> 3<sup>Lumen</sup> 16<sup>cm</sup><sub>catheter</sub>



**Arrowg+ard Blue Plus® Three-Lumen CVC**