





REF CA-25123-F

20 cm catheter length .035 inch dia. spring-wire quide

Arrowg+ard Blue® Large-Bore Three-Lumen CVC Kit for High Volume Infusions

Contents:

- Arrowg+ard Blue® Three-Lumen Catheter: 12 Fr. (4.1 mm OD) x 20 cm Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End "J" Tip on Other) with Arrow Advancer
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle

- Pressure Transduction Probe
 Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
 Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction
- Syringe Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe

- Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
 Dust Cap: Non-Vented
 SharpsAway® Disposal Cup
 Drape: 24" x 38" (61 cm x 91 cm) with 4" (10 cm) fenestration

- Drape: 24" x 30" (b) GIII x 91 GII) with 7 (10.5 Scalpel: #11 Gauze Pad: 2" x 2" (5 cm x 5 cm) Gauze Pad: 4" x 4" (10 cm x 10 cm) Suture: 3-0 Braided Silk with Straight Needle

| Rx only |
|--|
| Warning: Read all package insert warnings, precautions, and instructions |
| prior to use. Failure to do so may result in severe patient injury or death. |

California Prop. 65

MARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

| Lumen | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) | |
|-------------------|----------------------------|-------------------------------------|--|
| Distal (16 Ga.) | 0.4 | 2457 | |
| Medial (12 Ga.) | 1.3 | 10151 | |
| Proximal (12 Ga.) | 1.2 | 11039 | |

^{*} Priming volumes are approximate and are done without accessories.










Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL051425 R01 (2022-02)













 $[\]ensuremath{\dagger}$ Flow rate values are approximate and are determined using deionized water at 100 cm head height.