



REF CA-22702 (IPNXXXXXX)



7 Fr. 2 Lumen 16 cm catheter length .025 inch dia. spring-wire guide

# Arrowgard Blue® Two-Lumen CVC

**Contents:**

- 1: Arrowgard Blue® Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Pressure Transduction Probe origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 2: Dust Cap: Non-Vented origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: SecondSite™ Adjustable Hub: Fastener origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: SecondSite™ Adjustable Hub: Catheter Clamp origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® Disposal Cup origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Scalpel: #11 origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Suture: 3-0 Braided Silk with Straight Needle origin:xxxxxxxxxxxxxxxxxxxxxxxx

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**

**WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.**

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

**Contraindications: The Arrowgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.**



Lumen

Priming Volume\* (mL)

Gravity Flow Rate† (mL/hr)

Distal (18 Ga.)	0.33	1079
Proximal (14 Ga.)	0.50	4648

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)LotNumber

Set Product of xxx  
Packaged in xxx

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA



LBL051439 R01 (2024-05)



REF CA-22702

7 Fr. 2 Lumen

16 cm catheter length

.025 inch dia. spring-wire guide

# Arrowgard Blue® Two-Lumen CVC

