





Arrowg+ard Blue® Three-Lumen CVC

- Arrowg+ard Blue® Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
- Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One 1. End - "J" Tip on Other) with Arrow Advancer
- Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- Syringe: 5 mL Luer-Slip 1.
- Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- Dust Cap: Non-Vented 3:
- SecondSite™ Adjustable Hub: Fastener
- SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

cs da ni et
fi fr de el
hu is it lv
It no pl pt
ro ru sr sk
sl es sv tr
uk
$\underline{\mathring{\mathbb{N}}} \underline{\mathbb{M}} \underline{\mathbb{M}} \underline{\mathbb{M}}$
A 2 STERBUZE
STERILE EO *
25°C (77°F)

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (16 Ga.)	0.42	2656	
Medial (18 Gá.)	0.38	1353	
Proximal (18 Ga.)	0.41	1477	



[†] Flow rate values are approximate and are determined using deionized water at 100 cm head





EU Authorized Representative and Importer:

Teleflex Medical EC REP

IDA Business and Technology Park Dublin Road, Athlone 2797 Co. Westmeath, Ireland

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL051510 R02 (2023-11)

















