





Arrowg+ard Blue® Four-Lumen CV0

Contents:

Arrowg+ard Blue® Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 16 cm Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark

Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW

Syringe: 5 mL Luer-Slip

Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm

Dust Cap: Non-Vented
SecondSite™ Adjustable Hub: Fastener
SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C

(104°F). Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

de el hu	\bigcirc	
de el hu	4	is
it Iv It	X	no
pl pt ro		ru
sr sk sl		es
sv tr uk		
/ MD		$\widetilde{f i}$
A 2) (2 ERINZE
STERILE EO (\supset
类学		8
LAXXEX	(77 (77	°C °F)

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (16 Ga.)	0.40	2415	
Medial 1 (14 Ga.)	0.48	4723	
Medial 2 (18 Ga.)	0.37	1311	
Proximal (18 Ga.)	0.39	1403	



† Flow rate values are approximate and are determined using deionized water at 100 cm head height.





EU Authorized Representative and Importer:

EC REP

IDA Business and Technology Park Dublin Road, Athlone 2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

> Product of xxx Packaged in xxx



www.teleflex.com

(01)XXXXXXXXXXXXXXX

(17)YYMMDD (11)YYMMDD

(10)Lot Number

LBL051514 R03 (2025-01)





OT Lot Number

Y-MM-DD



Lumen

YY-MM-DD

catheter lenath

\$\$\ \text{1032} \ \text{spring-wire} \ \text{guide}





