





## Arrowg+ard Blue® Four-Lumen CV0

## Contents:

Arrowg+ard Blue® Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 20 cm Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End-"J" Tip on Other) with Arrow Advancer with ECG Mark Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW

Syringe: 5 mL Luer-Slip

Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm

Dust Cap: Non-Vented
SecondSite™ Adjustable Hub: Fastener
SecondSite™ Adjustable Hub: Catheter Clamp All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C

(104°F). Fluid path components are non-pyrogenic. Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

de el hu is
it ly lt no
pl pt ro ru
sr sk sl es
<b>—</b>
MD i
STERILE EO
<b>**                                   </b>
LAXEX (77°F)

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (16 Ga.)	0.45	2108	
Medial 1 (14 Ga.)	0.58	3860	
Medial 2 (18 Ga.)	0.39	1096	
Proximal (18 Ga.)	0.40	1129	



† Flow rate values are approximate and are determined using deionized water at 100 cm head height.





IDA Business and Technology Park Dublin Road, Athlone 2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

> Product of xxx Packaged in xxx




www.teleflex.com

(01)XXXXXXXXXXXXXXX

(17)YYMMDD (11)YYMMDD

(10)Lot Number

OT Lot Number







LBL051526 R03 (2024-10)





Lumen

catheter

.032 spring-wire guide





