





## Arrowg+ard Blue® Two-Lumen CV(

## Contents:

Arrowg+ard Blue® Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW

Pressure Transduction Probe

Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm

Dust Cap: Non-Vented
SecondSite™ Adjustable Hub: Fastener
SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU California Prop. 65

MARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.







Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (16 Ga.) Proximal (16 Ga.)	0.45 0.49	2540 2888	
* Priming volumes	are approximate	e and are done with	out accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.





**EU Authorized** Representative and Importer:

Teleflex Medical EC REP IDA Business and Technology Park Dublin Road, Athlone

2797 Co. Westmeath, Ireland All components are CE 2797 unless otherwise noted.

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA


> Product of xxx Packaged in xxx




(01)XXXXXXXXXXXXXXX (17)YYMMDD

(11)YYMMDD (10)Lot Number

LBL051528 R03 (2025-01)















Lumen



spring-wire





