





(IPNXXXXXX)

5 cm catheter length

.021 inch dia. spring-wire guide

## Small Vein Arrowg+ard Blue® Two-Lumen CVC

- <sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Injury or death. www.teletiex.com/iru
California Prop. 65

MARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov
Not made with natural rubber latex.
Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents


Fluid path components are non-pyrogenic

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (18 Ga.)	0.26	1760	
Proximal (20 Ga.)	0.30	1253	


\* Priming volumes are approximate and are done without accessories . ....... это ате арргохітнате and are done without accessories. † Flow rate values are approximate and are determined using deionized water at 100 cm head height.













(01)XXXXXXXXXXXXXXX (17)YYMMDD (11)YYMMDD



Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

**Teleflex**®

LBL052065 R04 (2025-02)







5 cm catheter length

.021 inch dia. spring-wire quide











