## PB-10000-850 8.5 X 10.5



1: Checklist/CLIP Sheet

1: Sterile Procedure Sign

Medication Label: 1% Lidocaine

Gauze Pad: 4" x 4" (10 cm x 10 cm)

Surgical Apparel: Impervious Gown

Suture: 3-0 Silk with Curved Needle

<sup>1</sup>A trademark of Becton, Dickinson and Company.

Surgical Apparel: Bouffant Cap Surgical Apparel: Mask with Eye Shield

1: HemoHopper® Fluid Receptacle

<sup>3</sup>A registered trademark of 3M Company.

<sup>5</sup>A registered trademark of C. R. Bard, Inc.

Dressing: Tegaderm®3 10 cm x 12 cm

Dressing: STATLOCK®5 CV Plus Catheter Stabilization Device

<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

<sup>4</sup>A registered trademark of Johnson & Johnson Corporation

2: Gauze Pad: 2" x 2" (5 cm x 5 cm)

Dressing: BIOPATCH®4

1:

5:

1:

1.

1:

1:

1:

1.

1:

## Contents:

- 1: Arrowg+ard Blue® Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 8 cm 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) PTFE-Coated
- (Straight Soft Tip on One End "J" Tip on Other) with Arrow Advancer
  Catheter: 20 Ga. x 1.16" (30 mm) Insyte™1 Autoguard™ Radiopaque over 22 Ga. Introducer Needle
- 1: Injection Needle: SafetyGlide<sup>™1</sup> 25 Ga. x 1" (2.54 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW and 5 mL Luer-Slip Syringe
- 1: Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 2: 10 mL Luer-Lock Pre-Filled Saline Syringe
- 1: 5 mL 1% Lidocaine HCI Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep<sup>®2</sup> One-Step Solution with Hi-Lite Orange<sup>™</sup> Tint
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite<sup>™</sup> Adjustable Hub: Catheter Clamp
- 1: SharpsAway® II Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Towel
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1. Salety Scalpel. #11

## Not made with natural rubber latex.

20°C (68°F)

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (20 Ga.)	0.3	1121	
Medial (22 Ga.)	0.4	587	
Proximal (22 Ga.)	0.4	896	

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

LBL052891 R02 (2024-05)