



One-Lumen CVC


1: One-Lumen Catheter: 16 Ga. (1.8 mm OD) x 20 cm origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Spring-Wire Guide: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Injection Needle: 25 Ga. x 1" (2.54 cm) origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Luer-Slip Syringe origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Injection Needle: 20 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Syringe: 3 mL Luer-Lock origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®1 One-Step Solution with Hi-Lite Orange™ Tint
origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: SecondSite™ Adjustable Hub: Fastener origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: SecondSite™ Adjustable Hub: Catheter Clamp origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: SharpsAway® Disposal Cup origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Drape: 26" x 18" (66 cm x 46 cm) with adhesive origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Scalpel: #11 origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Gauze Pad: 4" x 4" (10 cm x 10 cm) origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Suture: 3-0 Nylon with Straight Needle origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx
1: Tubing±: 6" Extension origin:xxxxxxxxxxxxxxxxxxxxxxxxxxxxx

1A registered trademark of CareFusion or one of its subsidiaries.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

3179

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



LBL053204 R01 (2024-05)

