## PB-10000-850 8.5 X 10.5



REF ASK-04300-SPH (IPNXXXXXX)

Lumen

Fr.

5





## **One-Lumen CVC Kit**

## Contents:

- 1: One-Lumen Catheter: 5 Fr. (1.8 mm OD) x 20 cm
- Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) 1: (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW 1. Introducer Needle
- Pressure Transduction Probe 1.
- 1. Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- Injection Needle: Eclipse<sup>™1</sup> 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe Injection Needle: Eclipse<sup>™1</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock 1:
- Syringe
- Syringe: 10 mL Luer-Lock 1.
- Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm 1:
- SecondSite™ Adjustable Hub: Fastener 1:
- SecondSite™ Adjustable Hub: Catheter Clamp 1:
- SharpsAway® II Locking Disposal Cup 1:
- SharpsAway® Disposal Cup
- Maximal Barrier Drape™ with 4" fenestration 1.
- Towel 1:
- Needle Holder 1:
- Safety Scalpel: #11 1.

LÀTEX

- Checklist/CLIP Sheet 1:
- Sterile Procedure Sign 1:
- 2: Gauze Pad: 3" x 3" (7.6 cm x 7.6 cm)

## Not made with natural rubber latex.

- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1:
- Surgical Apparel: Impervious Gown Dressing: I.V. Advanced 8.5 cm x 11.5 cm 1:
- Surgical Apparel: Bouffant Cap 1:
- Surgical Apparel: Mask with Eye Shield 1:
- Suture: 3-0 Silk with Curved Needle 1:
- 1: Valve: MaxZero<sup>™2</sup> Needleless Connector
- 1: HemoHopper® Fluid Receptacle
- <sup>1</sup>A trademark of Becton, Dickinson and Company.
- <sup>2</sup>A trademark of CareFusion or one of its subsidiaries

O	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (16 Ga.)	0.4	2498	

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

LBL053363 R02 (2024-04)



ASK ARRON