





1: Scalpel: #11
2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.

13 cm catheter length .021 inch dia. spring-wire guide

Small Vein Two-Lumen CVC

Contents:

- 1: Two-Lumen Catheter: 5 Fr. (1.8 mm OD) x 13 cm
 1: Spring-Wire Guide: .021" (0.53 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J" Tip on Other)
 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
 1: Introducer Needle: 25 Ga. x 1-1/2" (3.81 cm) XTW
 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) XTW
 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle with 5 mL Luer-Slip Syringe
 1: Syringe: 3 mL Luer-Lock
 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
 1: 5 mL 1% Lidocaine HCI Solution and Alcohol Prep
 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint
 2: Dust Cap: Non-Vented
 1: SecondSite™ Adjustable Hub: Fastener
 1: SecondSite™ Adjustable Hub: Catheter Clamp
 1: SharpsAway® Disposal Cup
 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration

| 6 | |
|-------|--|
| Lumon | |

Distal (18 Ga.) Proximal (20 Ga.)

* Priming volumes are approximate and are done without accessories.
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

Store between 20 - 25°C (68 - 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents

Fluid path components are non-pyrogenic.

























Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL053644 R02 (2022-06)





13 cm catheter length

.021 inch dia. spring-wire guide

