



REF AK-16553

(IPNXXXXXX)

## 5.5<sup>Fr.</sup>

### 3 Lumen

**13** cm  
catheter  
length

**.018** inch dia.  
spring-wire  
guide

## Three-Lumen CVC

## Contents:

- 1: Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 13 cm origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)  
origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Injection Needle: 25 Ga. x 1" (2.54 cm) origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe  
origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Syringe: 3 mL Luer-Lock origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm origin:xxxxxxxxxxxxxxxxxxxxxx  
1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep origin:xxxxxxxxxxxxxxxxxxxxxx  
1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® 1 One-Step Solution with Hi-Lite Orange™ Tint  
origin:xxxxxxxxxxxxxxxxxxxxxx  
3: Dust Cap: Non-Vented origin:xxxxxxxxxxxxxxxxxxxxxx  
1: SecondSite™ Adjustable Hub: Fastener origin:xxxxxxxxxxxxxxxxxxxxxx  
1: SecondSite™ Adjustable Hub: Catheter Clamp origin:xxxxxxxxxxxxxxxxxxxxxx  
1: SharpsAway® Disposal Cup origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Scalpel: #11 origin:xxxxxxxxxxxxxxxxxxxxxx  
2: Gauze Pad: 2" x 2" (5 cm x 5 cm) origin:xxxxxxxxxxxxxxxxxxxxxx  
5: Gauze Pad: 4" x 4" (10 cm x 10 cm) origin:xxxxxxxxxxxxxxxxxxxxxx  
1: Suture: 3-0 Braided Silk with Straight Needle origin:xxxxxxxxxxxxxxxxxxxxxx

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only**

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

**WARNING:** Cancer and Reproductive Harm.  
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



Lumen

Priming  
Volume\*  
(mL)Gravity  
Flow  
Rate†  
(mL/hr)

Distal (20 Ga.)	0.32	988
Medial (22 Ga.)	0.38	370
Proximal (22 Ga.)	0.34	422

\* Priming volumes are approximate and are done without accessories.


† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



**LOT** Lot Number



(01)XXXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

 Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Set Product of xxx  
Packaged in xxx

**Teleflex®**  
www.teleflex.com

LBL053648 R02 (2024-04)



**REF** **AK-16553**

## 5.5<sup>Fr.</sup>

### 3 Lumen

**13** cm catheter length

**.018** inch dia.  
spring-wire  
guide



## Three-Lumen CVC