



REF EU-47702-MSB (IPNXXXXXX)

ARROW

7 Fr. 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Arrowgard Blue Plus® Two-Lumen CVC

Contents:

- 1: Arrowgard Blue Plus®1 Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Injection Needle: Eclipse™2 23 Ga. x 1" (2.54 cm)
- 1: Injection Needle: Eclipse™2 25 Ga. x 1" (2.54 cm) RW
- 1: Pressure Transduction Probe
- 1: Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
- 1: Syringe: 5 mL Luer-Lock
- 1: Syringe: 3 mL Luer-Lock
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® II Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Maximal Barrier Drape: 55" x 94-1/2" (140 cm x 240 cm) with 4" (10 cm) fenestration, with adhesive and covered perforation
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 2: Stopcock: 4-way High-Flow
- 1: Sterile Procedure Sign
- 2: Prep Sponge Swab
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 6: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Surgical Apparel: Cap
- 1: Surgical Apparel: Impervious Gown
- 1: Surgical Apparel: Mask with Shield
- 1: Suture: (0) Silk with Curved Needle
- 1: Tray: Prep

1Licensed under US Patent No. 7,329,412.
2A trademark of Becton, Dickinson and Company.
All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).
Fluid path components are non-pyrogenic.

Contraindications: The Arrowgard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.




Lumen

Priming Volume* (mL)

Gravity Flow Rate† (mL/hr)

Distal (14 Ga.)	0.57	4731
Proximal (18 Ga.)	0.40	913

* Priming volumes are approximate and are done without accessories.
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



LOT Lot Number
YYYY-MM-DD
YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

EU Authorized Representative and Importer: Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LBL055046 R03 (2025-03)

