







20 cm catheter length .018 inch dia. spring-wire quide

Pressure Injectable One-Lumen Midline

Contents:

- TaperFree® Catheter: 4 Fr.(1.45 mm OD) x 20 cm Pressure Injectable Radiopaque Polyurethane with T-Port Connector, Extension Line Clamp, Contamination Guard and Stiffening Stylet GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dileter
- Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
- Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
- Syringe: 10 mL Luer-Lock SharpsAway® II Locking Disposal Cup Catheter Trimmer
- Safety Scalpel: #11
- Paper Tape Measure
- Tourniquet
 Dressing: GRIP-LOK®1 Catheter Fixation Device
- Valve: MicroCLAVE®2 Neutral Displacement Connector
- ¹A registered trademark of TIDI.
- ²A registered trademark of ICU Medical, Inc.

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	Pump Flow Rate†† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Yellow (18 Ga.)	0.37	1880	9324	5

^{*} Priming volumes are approximate and are done without accessories.

^{††} Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.
** Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 152 cm administration set tubing.





LOT













Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL055050 R02 (2022-01)















 $[\]dagger$ Flow rate values are approximate and are determined using deionized water at 100 cm head height.