
ARR(





(IPNXXXXXX)

8 cm catheter length

Two-Lumen CVC

Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 8 cm Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer

Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle

Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW

Syringe: 5 mL Luer-Slip Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm

Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
SecondSite™ Adjustable Hub: Catheter Clamp
SecondSite™ Adjustable Hub: Fastener

SharpsAway® II Locking Disposal Cup

Safety Scalpel: #11

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex.

Fluid path components are non-pyrogenic.


















































Distal (22 Ga.) Proximal (22 Ga.) 0.28 548 880

* Priming volumes are approximate and are done without accessories. † Flow rate values are approximate and are determined using deionized water at 100 cm head height.



EU Authorized Representative and Importer:



Teleflex Medical IDA Business and Technology Park Dublin Road, Athlone 2797 Co. Westmeath, Ireland


All components are CE 2797 unless otherwise noted.

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

Product of xxx Packaged in xxx

LOT LotNumber







(01)XXXXXXXXXXXXXXX (17)YYMMDD (11)YYMMDD (10)LotNumber

)feleflex® www.teleflex.com

LBL055286 R02 (2024-08)



.018 inch dia. spring-wire guide



