**ARRC** 

Product of XXX Packaged in XXX





(IPNXXXXXX)

3 cm catheter length

■018 inch dia. spring-wire guide

## **Pediatric Two-Lumen CVC**

- Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 13 cm Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J" Tip on Other) with Arrow Advancer Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- Syringe: 5 mL Luer-Slip
- Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm

- Dust Cap: Non-Vented
  SecondSite™ Adjustable Hub: Catheter Clamp
- SecondSite™ Adjustable Hub: Fastener

All components are CE 2797 unless otherwise noted

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

▲ WARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex. Fluid path components are non-pyrogenic.



















Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (22 Ga.)	0.27	363	
Proximal (22 Ga.)	0.31	277	

Priming volumes are approximate and are done without accessories.



Teleflex Medical EC REP IDA Business and E 2797 Technology Park Dublin Road, Athlone Co. Westmeath, Ireland


Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA




LBL055871 R01 (2024-02)

**ARROW** 









(01)XXXXXXXXXXXXXXX (11)YYMMDD (10)LotNumber



**.**018 inch dia. spring-wire guide





<sup>†</sup> Flow rate values are approximate and are determined using deionized water at 100 cm head height.