





20 cm catheter length

.032 inch dia. spring-wire guide

Pressure Injectable Three-Lumen CVC

- Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm, Pressure Iniectable
- Injectation
 Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x
 23-5/8" (60 cm) (Straight Soft Tip on One End "J" Tip on Other) with
 Arrow Advancer with ECG Mark
 Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW
 Introducer Needle
- Introducer Needle
 Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
 Pressure Transduction Probe
 Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
 Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
 Dust Cap: Non-Vented
 SecondSite™ Adjustable Hub: Fastener
 SecondSite™ Adjustable Hub: Catheter Clamp

- Dressing: GRIP-LOK®1 Catheter Fixation Device
- ¹A registered trademark of TIDI.

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (16 Ga.)	0.47	2656	10
Medial (18 Ga.)	0.40	1245	5
Proximal (18 Ga.)	0.42	1328	5

^{*} Priming volumes are approximate and are done without accessories



LOT



















Teleflex Medical IDA Business and Technology Park
Dublin Road, Athlone
7797 Co. Westmeath, Ireland

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL056286 R01 (2023-09)









1032 inch dia. spring-wire guide







[†] Flow rate values are approximate and are determined using deionized water at 100 cm head height.

** Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 152 cm administration set tubing.