



REF PK-01451-WNM (IPNXXXXXX)

4 Fr. 1 Lumen 55 cm catheter length

# One-Lumen PICC

**Contents:**

- 1: One Lumen Catheter: 4 Fr. (1.45 mm OD) x 55 cm, T-Port Connector, Blue FlexTip®, Contamination Guard and Placement Wire
- 1: Peel-Away Sheath: 15 Ga. x 1-1/4" (3.18 cm) Radiopaque over 17 Ga. TW Introducer Needle
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 10 mL Luer-Slip
- 1: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SharpsAway® Disposal Cup
- 1: Towel
- 1: Towel: 24" x 36" (61 cm x 91 cm)
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Patient ID Card
- 1: Chart Sticker
- 1: Patient Information Booklet
- 2: Paper Tape Measure
- 1: Tourniquet
- 1: Dressing: OpSite®1 14 cm x 10 cm
- 5: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 6: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: STATLOCK®2 PICC Plus Catheter Stabilization Device
- 2: Surgical Apparel: Mask
- 1: Tape: Steri-Strip®3

<sup>3</sup>A registered trademark of 3M Company.

<sup>1</sup>A registered trademark of T.J. Smith & Nephew, Limited.  
<sup>2</sup>A registered trademark of C. R. Bard, Inc.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**

**Not made with natural rubber latex.**

Fluid path components are non-pyrogenic.

**Contraindications: The PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.**



Lumen

Priming Volume\* (mL)

Gravity Flow Rate† (mL/hr)

Pump Flow Rate‡ (mL/hr)

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga.)	0.49	568	4405

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.



LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)LotNumber

Product of xxx  
Packaged in xxx

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA



LBL056629 R02 (2024-07)



REF PK-01451-WNM

4 Fr. 1 Lumen 55 cm catheter length

# One-Lumen PICC

