

Dimensions: 6-3/8" (16 cm) x 8-1/4" (21 cm)

ARROW®

REF **PI-01351-SS** (IPNXXXXXX)for use
with
Arrow

3

Fr. PICC

3

cm

.016

inch dia.
spring-wire
guide

Insertion Set for use with 3 Fr. PICC

Contents:

- 1: GlideThru™ Peel-Away Sheath: 3 Fr. x 1-1/4" (3 cm) Radiopaque over 3 Fr. Dilator
origin:xxxxxxxxxxxxxxxxxxxxxx
- 1: Spring Wire Guide Assembly: .016" (0.40 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip On
Both Ends) origin:xxxxxxxxxxxxxxxxxxxxxx
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
origin:xxxxxxxxxxxxxxxxxxxxxx
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW origin:xxxxxxxxxxxxxxxxxxxxxx

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex. Fluid path components are non-pyrogenic.



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Set Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LOT Lot Number

YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL056668 R01 (2024-05)