



REF PI-01351-SS

for use
with
Arrow

3

Fr. PICC

3

cm

.016 inch dia.
spring-wire
guide

Insertion Set for use with 3 Fr. PICC

Contents:

- 1: GlideThru™ Peel-Away Sheath: 3 Fr. x 1-1/4" (3 cm) Radiopaque over 3 Fr. Dilator
- 1: Spring Wire Guide Assembly: .016" (0.40 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip On Both Ends)
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use.


Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



LOT



 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL056668 R00 (2022-04)