



REF CDC-16123-XP1A

12 Fr. 3 Lumen 25 cm catheter length .035 inch dia. spring-wire guide

ARROW

ErgoPack

## Pressure Injectable Three-Lumen Hemodialysis Catheterization for High Volume Infusions

### Contents:

- |  |  |
|--|--|
| <ul style="list-style-type: none"> <li>1: Three-Lumen Catheter: 12 Fr. (4.06 mm OD) x 25 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps</li> <li>1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer</li> <li>1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle</li> <li>1: Injection Needle: SafetyGlide™ 1 25 Ga. x 1" (2.54 cm)</li> <li>1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)</li> <li>1: Pressure Transduction Probe</li> <li>1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe</li> <li>1: Injection Needle: SafetyGlide™ 1 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe</li> <li>1: Syringe: 3 mL Luer-Lock</li> <li>1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm</li> <li>1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm</li> <li>2: 10 mL Luer-Lock Pre-Filled Saline Syringe</li> <li>1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep</li> <li>1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® 2 One-Step Solution with Hi-Lite Orange™ Tint</li> <li>1: SharpsAway® Disposal Cup</li> </ul> | <ul style="list-style-type: none"> <li>1: SharpsAway® II Locking Disposal Cup</li> <li>1: Maximal Barrier Drape™ with 4" fenestration</li> <li>1: Towel</li> <li>1: Needle Holder</li> <li>1: Safety Scalpel: #11</li> <li>1: Checklist/CLIP Sheet</li> <li>1: Sterile Procedure Sign</li> <li>1: Medication Label: 1% Lidocaine</li> <li>2: Gauze Pad: 2" x 2" (5 cm x 5 cm)</li> <li>5: Gauze Pad: 4" x 4" (10 cm x 10 cm)</li> <li>1: Dressing: Tegaderm® 10 cm x 12 cm</li> <li>1: Surgical Apparel: Impervious Gown</li> <li>1: Dressing: BIOPATCH® 4</li> <li>1: Surgical Apparel: Bouffant Cap</li> <li>1: Surgical Apparel: Mask with Eye Shield</li> <li>1: Suture: 3-0 Silk with Curved Needle</li> <li>1: HemoHopper® Fluid Receptacle</li> </ul> <p><sup>1</sup>A trademark of Becton, Dickinson and Company.<br/> <sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.<br/> <sup>3</sup>A registered trademark of 3M Company.<br/> <sup>4</sup>A registered trademark of Johnson &amp; Johnson Corporation.</p> |
|--|--|

### Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

**Contraindications: The Multi-Lumen Acute Hemodialysis Catheter (AHDC) is not designed for long-term (≥ 30 days) hemodialysis or for insertion into thrombosed vessels.**

