



REF CDC-25142-1A

14 Fr. 2 Lumen 20 cm catheter length .035 inch dia. spring-wire guide

ARROW

ErgoPack®

# Arrowg+ard Blue® Two-Lumen Hemodialysis Catheterization for High Volume Infusions


**Contents:**

- 1: Two-Lumen Catheter: 14 Fr. (4.72 mm OD) x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment, Extension Line Clamps
  - 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
  - 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
  - 1: Injection Needle: SafetyGlide™ 1 25 Ga. x 1" (2.54 cm)
  - 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
  - 1: Pressure Transduction Probe
  - 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
  - 1: Injection Needle: SafetyGlide™ 1 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
  - 1: Syringe: 3 mL Luer-Lock
  - 1: Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm
  - 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
  - 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
  - 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Solution with Hi-Lite Orange™ Tint
  - 1: SharpsAway® Disposal Cup
  - 1: SharpsAway® II Locking Disposal Cup
  - 1: Maximal Barrier Drape™ with 4" fenestration
  - 1: Needle Holder
  - 1: Safety Scalpel: #11
  - 1: Checklist/CLIP Sheet
  - 1: Sterile Procedure Sign
  - 1: Medication Label: 1% Lidocaine
  - 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
  - 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
  - 1: Surgical Apparel: Impervious Gown
  - 1: Surgical Apparel: Bouffant Cap
  - 1: Surgical Apparel: Mask with Eye Shield
  - 1: Suture: 3-0 Braided Silk with Straight Needle
  - 1: Suture: 3-0 Silk with Curved Needle
  - 1: HemoHopper® Fluid Receptacle
- <sup>1</sup>A trademark of Becton, Dickinson and Company.  
<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Not made with natural rubber latex.**


Store between 20 – 25°C (68 – 77°F). Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

**Contraindications:** The Arrow large-bore two-lumen catheter is not designed for long-term (≥ 30 days) hemodialysis or for use in patients with thrombosed vessels. The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs.



Priming Volume*	
Distal	1.70 mL
Proximal	1.50 mL

\*Priming volumes are approximate and are done without accessories.



† In vitro tests performed using 3:1 glycerine/saline at 37°C.

