





REF AK-15123-PCDC

20 cm catheter length

.035 inch dia. spring-wire guide

Pressure Injectable Three-Lumen Hemodialysis Catheterization for High Volume **Infusions**

Contents:

- Three-Lumen Catheter: 12 Fr. (4.06 mm OD) x 20 cm Pressure
- Inrec-Lumen Catheter: 12 Fr. (4.06 mm OD) x 20 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps

 Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End "J" Tip on Other) with Arrow
- Advancer Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW
- Introducer Needle
 Injection Needle: 25 Ga. x 1" (2.54 cm)
 Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

- Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
 Pressure Transduction Probe
 Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow®
 Raulerson Spring-Wire Introduction Syringe
 Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip
- Syringe Syringe: 3 mL Luer-Lock
- Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®1 One-Step
- Solution with Hi-Lite Orange™ Tint SharpsAway® II Locking Disposal Cup

- SharpsAway® Disposal Cup Drape: 41" x 55" with 4" fenestration Needle Holder Safety Scalpel: #11

- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
 1: Suture: 3-0 Silk with Curved Needle
- ¹A registered trademark of CareFusion or one of its subsidiaries.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU California Prop. 65

____ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov
Not made with natural rubber latex.

Store between 20 - 25°C (68 - 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Contraindications: The Multi-Lumen Acute Hemodialysis Catheter (AHDC) is not designed for long-term (≥ 30 days) hemodialysis or for insertion into thrombosed vessels.







Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL057102 R01 (2023-08)









