



REF AK-12123-PCDC

12 Fr. 3 Lumen 16 cm catheter length .035 inch dia. spring-wire guide

Pressure Injectable Three-Lumen Hemodialysis Catheterization for High Volume Infusions

Contents:

- | | |
|--|---|
| 1: Three-Lumen Catheter: 12 Fr. (4.06 mm OD) x 16 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps | 1: SharpsAway® Disposal Cup |
| 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: Drape: 41" x 55" with 4" fenestration |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: Needle Holder |
| 1: Injection Needle: SafetyGlide™ 1 25 Ga. x 1" (2.54 cm) | 1: Safety Scalpel: #11 |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Pressure Transduction Probe | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: Suture: 3-0 Silk with Curved Needle |
| 1: Injection Needle: SafetyGlide™ 1 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | ¹ A trademark of Becton, Dickinson and Company. |
| 1: Syringe: 3 mL Luer-Lock | ² A registered trademark of CareFusion or one of its subsidiaries. |
| 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm | |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | |
| 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep® One-Step Solution with Hi-Lite Orange™ Tint | |
| 1: SharpsAway® II Locking Disposal Cup | |

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

WARNING: Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Contraindications: The Multi-Lumen Acute Hemodialysis Catheter (AHDC) is not designed for long-term (≥ 30 days) hemodialysis or for insertion into thrombosed vessels.



LOT



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL057104 R01 (2023-08)



REF AK-12123-PCDC

12 Fr. 3 Lumen 16 cm catheter length .035 inch dia. spring-wire guide

Pressure Injectable Three-Lumen Hemodialysis Catheterization for High Volume Infusions