





## Two-Lumen CVC for use only with Arrow MAC™ Two-Lumen Central Venous Access Device

## Contents:

Two-Lumen Catheter: 7 Fr. (2.5 mm OD) Catheter Contamination Guard

Dust Cap: Non-Vented

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.







6 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (16 Ga.)	0.37	2664	
Proximal (16 Ga.)	0.45	2913	

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.







All components are CE 2797 unless otherwise noted.

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

> Product of xxx Packaged in xxx

1feleflex® www.teleflex.com

(01)XXXXXXXXXXXXXXX (17)YYMMDD

(11)YYMMDD (10)Lot Number







LBL057194 R03 (2024-09)



Lumen

cm catheter





