



REF ML-00702



Two-Lumen CVC for use only with Arrow MAC™ Two-Lumen Central Venous Access Device

Contents: _____

- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) Catheter Contamination Guard
- 2: Dust Cap: Non-Vented

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.37	2664
Proximal (16 Ga.)	0.45	2913

* Priming volumes are approximate and are done without accessories.
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

CAS # 7440-48-4

EU Authorized Representative and Importer:

EC REP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
 2797

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL057194 R02 (2023-10)

LOT



7 Fr. 2 Lumen 16 cm catheter length with MAC

REF ML-00702

Two-Lumen CVC for use only with Arrow MAC™ Two-Lumen Central Venous Access Device