



REF CS-15142-F

20 cm catheter length .038 inch dia. spring-wire guide

Two-Lumen Hemodialysis Catheter for High Volume Infusions

Contents:

- Two-Lumen Catheter: 14 Fr. (4.72 mm OD) x 20 cm Radiopaque Polyurethane with Blue 1: FlexTip®, Extension Line Clamps
- FIEXTIP®, Extension Line Clamps
 Spring-Wire Guide, Marked: .038" (0.97 mm) dia. x 27-1/2" (70 cm) (Straight Soft Tip on One End -"J" Tip on Other) with Arrow Advancer
 Catheter: 16 Ga. x 1-7/8" (4.8 cm) Radiopaque over 19 Ga. RW Introducer Needle
 Pressure Transduction Probe
 Introducer Needle: 16 Ga. x 2.4 (2016) 6.5 Cm.)

- Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-
- Wire Introduction Syringe
 Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm
 Tissue Dilator: 16 Fr. (5.3 mm) x 11.4 cm
- Dust Cap: Non-Vented
 Drape: 24" x 36" (61 cm x 91 cm)
 Scalpel: #11

All components are CE 2797 unless otherwise noted.

Rx		

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU California Prop. 65

M WARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov
Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

Contraindications: The large-bore two-lumen catheter is not designed for long-term (≥ 30 days) hemodialysis or for use in patients with thrombosed vessels.

Lumen	Priming Volume* (mL)	Flow Rate (mL/min)	Venous Pressure Setting†† (mmHg)
Distal (10 Ga.)	1.7	250	100
Proximal (10 Ga.)	1.5	400 250 400	170 120 190

^{*} Priming volumes are approximate and are done without accessories













EU Authorized Representative and Importer:



Teleflex Medical IDA Business and Technology Park Dublin Road, Athlone
Co. Westmeath, Ireland

Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL057847 R01 (2023-08)











^{††} Venous Pressure settings are approximate and were determined by varying the blood flow rates from 250 to 400 mL/min in an in vitro experiment using blood analog solution.