



REF CS-25142-F

14<sup>Fr.</sup> 2<sup>Lumen</sup> 20<sup>cm</sup> catheter length .038<sup>inch dia.</sup> spring-wire guide

# Arrowg+ard Blue<sup>®</sup> Two-Lumen Hemodialysis Catheter for High Volume Infusions

**Contents:**

- 1: Two-Lumen Catheter: 14 Fr. (4.72 mm OD) x 20 cm Radiopaque Polyurethane with Blue FlexTip<sup>®</sup>, ARROWg+ard<sup>®</sup> Antimicrobial Surface Treatment, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .038" (0.97 mm) dia. x 27-1/2" (70 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 16 Ga. x 1-7/8" (4.8 cm) Radiopaque over 19 Ga. TW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow<sup>®</sup> Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm
- 1: Tissue Dilator 16 Fr. (5.3 mm) x 4 1/2" (11.4 cm)
- 2: Dust Cap: Non-Vented
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Scalpel: #11

All components are CE 2797 unless otherwise noted.

**Rx only**

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
California Prop. 65

**⚠ WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.**

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

**Contraindications:** The Arrow large-bore two-lumen catheter is not designed for long-term (≥ 30 days) hemodialysis or for use in patients with thrombosed vessels. The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs.



Lumen	Priming Volume* (mL)	Flow Rate (mL/min)	Venous Pressure Setting†† (mmHg)
Distal (10 Ga.)	1.7	250 400	100 170
Proximal (10 Ga.)	1.5	250 400	120 190

\* Priming volumes are approximate and are done without accessories.  
†† Venous Pressure settings are approximate and were determined by varying the blood flow rates from 250 to 400 mL/min in an in vitro experiment using blood analog solution.



LOT



EU Authorized Representative and Importer:

EC REP  
CE 2797  
Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex<sup>®</sup>**  
[www.teleflex.com](http://www.teleflex.com)

LBL057996 R00 (2021-06)



# Arrowg+ard Blue<sup>®</sup> Two-Lumen Hemodialysis Catheter for High Volume Infusions

REF CS-25142-F

14<sup>Fr.</sup> 2<sup>Lumen</sup> 20<sup>cm</sup> catheter length .038<sup>inch dia.</sup> spring-wire guide