



REF CS-26142-F

14^{Fr.} 2^{Lumen} 25^{cm} catheter length .038^{inch dia.} spring-wire guide

Arrowgard Blue® Two-Lumen Hemodialysis Catheter for High Volume Infusions

Contents:

- 1: Two-Lumen Catheter: 14 Fr. (4.72 mm OD) x 25 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowgard® Antimicrobial Surface Treatment, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .038" (0.97 mm) dia. x 27-1/2" (70 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 16 Ga. x 1-7/8" (4.8 cm) Radiopaque over 19 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm
- 1: Tissue Dilator 16 Fr. (5.3 mm) x 4 1/2" (11.4 cm)
- 2: Dust Cap: Non-Vented
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Scalpel: #11

All components are CE 2797 unless otherwise noted.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrow large-bore two-lumen catheter is not designed for long-term (≥ 30 days) hemodialysis or for use in patients with thrombosed vessels. The Arrowgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs.



Lumen	Priming Volume* (mL)	Flow Rate (mL/min)	Venous Pressure Setting†† (mmHg)
Distal (10 Ga.)	1.9	250	110
		400	180
Proximal (10 Ga.)	1.5	250	140
		400	200

* Priming volumes are approximate and are done without accessories.
 †† Venous Pressure settings are approximate and were determined by varying the blood flow rates from 250 to 400 mL/min in an in vitro experiment using blood analog solution.



EU Authorized Representative and Importer:
 Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
 CE 2797
 Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA
Teleflex®
www.teleflex.com
 LBL058001 R02 (2023-10)



Arrowgard Blue® Two-Lumen Hemodialysis Catheter for High Volume Infusions

REF CS-26142-F

14^{Fr.} 2^{Lumen} 25^{cm} catheter length .038^{inch dia.} spring-wire guide