



9 Fr. **8.5** cm sheath length **.035** inch dia. spring-wire guide

ARROW

PSI for use with 7 - 8 Fr. Catheters

Contents:

- | | |
|--|--|
| 1: Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 8.5 cm and Tissue Dilator | 1: Stopcock: 4-way High-Flow |
| 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) | 2: Prep Sponge Swab |
| (Straight Soft Tip on One End - "J" Tip on Other) with Arrow | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| Advancer | 1: Suture: 3-0 Braided Silk with Straight Needle |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW | 1: Tray: Prep |
| Introducer Needle | |
| 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) | |
| 1: Injection Needle: 25 Ga. x 1" (2.54 cm) | |
| 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW | |
| 1: Syringe: 3 mL Luer-Lock | |
| 1: Syringe: 5 mL Luer-Slip | |
| 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL | |
| 1: Obturator: 8 Fr. | |
| 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ | |
| Adapter | |
| 2: Dust Cap: Non-Vented | |
| 1: SharpsAway® Disposal Cup | |
| 1: Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with | |
| adhesive | |
| 1: Scalpel: #11 | |

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.


Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT



 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL058145 R00 (2021-12)



REF MTO-09903-TS

9 Fr. **8.5** cm sheath length **.035** inch dia. spring-wire guide

PSI for use with 7 - 8 Fr. Catheters

ARROW