C12122102A2 4.25' X 17.5'

DO NOT BEND





(IPNXXXXX)

20^{Ga.}

|REF| SAC-00520-PBX



Arterial Catheterization Set

Contents:

Catheter: 20 Ga. x 5 cm Radiopaque PEBA Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends) 1:


Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW All components are CE 2797 unless otherwise noted.


Catheter Priming Volume: 0.17 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Marking California Prop. 65 ▲ WARNING: Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex. Fluid path components are non-pyrogenic.





Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

Product of xxx Packaged in xxx



(10)Lot Number



LBL058185 R01 (2024-11)





Arterial Catheterization Set

ARROW DO NOT BEND