

DO NOT BEND

ARROW®



REF SAC-00522-PBX

(IPNXXXXXX)

22<sup>Ga.</sup>5<sup>cm</sup>  
catheter  
length

## Arterial Catheterization Set

## Contents:

- 1: Catheter: 22 Ga. x 5 cm Radiopaque PEBA  
 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)  
 1: Introducer Needle: 22 Ga. x 1-1/2" (4 cm) XTW

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.16 mL Priming volume is approximate.


Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.  
[www.teleflex.com/IFU](http://www.teleflex.com/IFU)

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

EU Authorized  
Representative and  
Importer:

EC REP



Teleflex Medical  
 IDA Business and  
 Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

Product of xxx  
Packaged in xxx

(01)XXXXXXXXXXXXXX  
 (17)YYMMDD  
 (11)YYMMDD  
 (10)Lot Number

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL058191 R02 (2025-03)



REF SAC-00522-PBX

22<sup>Ga.</sup>5<sup>cm</sup>  
catheter  
length

## Arterial Catheterization Set

ARROW®

DO NOT BEND