## DO NOT BEND







|REF| SAC-00818-PBX

(IPNXXXXXX)

18 Ga.

8 cm catheter length

## Arterial Catheterization Set

## Contents:

- Catheter: 18 Ga. x 8 cm Radiopaque PEBA Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- Introducer Needle: 18 Ga. x 2" (5 cm) XTW

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.19 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

California Prop. 65


Marning: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov Not made with natural rubber latex. Fluid path components are non-pyrogenic.



















(10)Lot Number



Teleflex Medical IDA Business and Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL058203 R01 (2024-01)





8 cm catheter length

## Arterial Catheterization Set

**ARROW** 

**DO NOT BEND**