## DO NOT BEND







|REF| SAC-01218

(IPNXXXXXX)

18<sup>Ga.</sup>

## Arterial Catheterization Set

## Contents:

- Catheter: 18 Ga. x 12 cm Radiopaque Polyurethane Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- Introducer Needle: 18 Ga. x 2" (5 cm) XTW

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.21 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU California Prop. 65

A WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex. Fluid path components are non-pyrogenic.

























All components are CE 2797 unless otherwise noted.



Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



(10)Lot Number

**Neleflex**® LBL058218 R02 (2025-03)



REF SAC-01218

18 Ga.

## Arterial Catheterization Set

**ARROW** 

**DO NOT BEND**