

DO NOT BEND



REF SAC-01618-PBX

(IPNXXXXXX)

18 Ga. 16 cm catheter length

# Arterial Catheterization Set

Contents:

- 1: Catheter: 18 Ga. x 16 cm Radiopaque PEBA
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.24 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



- en bg cs da nl et fi fr de el hu
- it lv lt no pl pt ro ru sk sl es
- sv tr

Product of xxx  
Packaged in xxx

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

EU Authorized Representative and Importer:

EC REP Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Teleflex®  
www.teleflex.com

LBL058239 R01 (2024-01)



REF SAC-01618-PBX

18 Ga. 16 cm catheter length

# Arterial Catheterization Set



DO NOT BEND