

DO NOT BEND

ARROW®



REF SAC-01620-PBX

(IPNXXXXXX)

20^{Ga.} 16^{cm catheter length}

Arterial Catheterization Set

Contents:

- 1: Catheter: 20 Ga. x 16 cm Radiopaque PEBA
 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 19-3/4" (50 cm) (Straight Soft Tip on Both Ends) with Arrow Advancer
 1: Introducer Needle: 20 Ga. x 2-3/4" (7 cm) XTW

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.21 mL Priming volume is approximate.


Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
www.teleflex.com/IFU

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



EU Authorized
Representative and
Importer:

EC REP



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

Product of xxx
Packaged in xxx



(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

Teleflex®
www.teleflex.com

LBL058245 R02 (2025-03)



REF SAC-01620-PBX

20^{Ga.} 16^{cm catheter length}

Arterial Catheterization Set

ARROW®

DO NOT BEND