

DO NOT BEND

ARROW®



REF SAC-02318-PBX
(IPNXXXXXX)

18^{Ga.} 23^{cm catheter length}

Arterial Catheterization Set

- Contents:**
- 1: Catheter: 18 Ga. x 23 cm Radiopaque PEBA
 - 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
 - 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW
- All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.28 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



en bg cs da nl et fi fr de el hu
it lv lt no pl pt ro ru sk sl es
sv tr

Product of xxx
Packaged in xxx

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number



EU Authorized
Representative and
Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL058251 R01 (2024-01)



REF SAC-02318-PBX

18^{Ga.} 23^{cm catheter length}

Arterial Catheterization Set

ARROW®

DO NOT BEND